






# ENDOCHOICE<sup>®</sup>

BONASTENT<sup>®</sup>


[www.endochoice.com](http://www.endochoice.com)

## Product Catalog 2010

Esophageal & Shim Anti-Reflux Valve  
Biliary  
M-Hilar  
Pyloric / Duodenal  
Colorectal  
Tracheal / Bronchial



# For The Quality of Living!




Polygon Mesh  
Surface Structure




Adapting Flexibility



Segmental  
Compression



Optimal Radial  
& Returning Force



Repositionable

## The Difference?



Hook Wire Structure





Cross Wire Structure




Hook & Cross Wire Structure (Polygon Mesh Surface)

- Low Rate of Shortening for Accurate Positioning
- Small Cavities for In-Growth Reduction
- Optimal Radial & Returning Force for Conformability
- Segmental Compression for Migration Reduction
- Reduced Delivery Device Diameter





# E Sophnaqaeal



Silicone Covering  
Retrieval Lasso  
Flange Ends for Migration Reduction  
Segmental Compression  
Repositionable

## Esophageal



Coding : B-BONASTENT / E-Esophageal / R-Reduced Shortening

| Order Number | Stent         |                    |                   | Delivery Device |                    |
|--------------|---------------|--------------------|-------------------|-----------------|--------------------|
|              | Diameter (mm) | Usable Length (mm) | Total Length (mm) | Diameter (mm)   | Usable Length (mm) |
| BE-1808      | 18            | 40                 | 80                | 5               | 700                |
| BE-1810      | 18            | 60                 | 100               | 5               | 700                |
| BE-1812      | 18            | 80                 | 120               | 5               | 700                |
| BE-1814      | 18            | 100                | 140               | 5               | 700                |
| BER-1808     | 18            | 40                 | 80                | 6               | 700                |
| BER-1810     | 18            | 60                 | 100               | 6               | 700                |
| BER-1812     | 18            | 80                 | 120               | 6               | 700                |
| BER-1814     | 18            | 100                | 140               | 6               | 700                |

## Esophageal Shim Anti-Reflux Valve




Coding : B-BONASTENT / E-Esophageal / V-Anti-Reflux Valve

| Order Number | Stent         |                    |                   | Delivery Device |                    |
|--------------|---------------|--------------------|-------------------|-----------------|--------------------|
|              | Diameter (mm) | Usable Length (mm) | Total Length (mm) | Diameter (mm)   | Usable Length (mm) |
| BEV-2209     | 22            | 50                 | 90                | 7               | 700                |
| BEV-2211     | 22            | 70                 | 110               | 7               | 700                |
| BEV-2213     | 22            | 90                 | 130               | 7               | 700                |
| BEV-2215     | 22            | 110                | 150               | 7               | 700                |

# Biliary





- Conformability
- Optimal Radial & Returning Force
- Non-Irritating Ends
- Segmental Compression
- Less Shortening
- Repositionable
- Silicone Covering

## Biliary



Coding : B-BONASTENT / B-Biliary

Delivery Device 1800mm Length for T.T.S.

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BB-100606    | 10            | 60          | 7               | 600                |
| BB-100706    | 10            | 70          | 7               | 600                |
| BB-100806    | 10            | 80          | 7               | 600                |
| BB-100906    | 10            | 90          | 7               | 600                |
| BB-100618    | 10            | 60          | 7               | 1800               |
| BB-100718    | 10            | 70          | 7               | 1800               |
| BB-100818    | 10            | 80          | 7               | 1800               |
| BB-100918    | 10            | 90          | 7               | 1800               |

## Biliary




Coding : B-BONASTENT / B-Biliary / C-Covered

Delivery Device 1800mm Length for T.T.S.

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BBC-100606   | 10            | 60          | 8               | 600                |
| BBC-100706   | 10            | 70          | 8               | 600                |
| BBC-100806   | 10            | 80          | 8               | 600                |
| BBC-100906   | 10            | 90          | 8               | 600                |
| BBC-100618   | 10            | 60          | 8               | 1800               |
| BBC-100718   | 10            | 70          | 8               | 1800               |
| BBC-100818   | 10            | 80          | 8               | 1800               |
| BBC-100918   | 10            | 90          | 8               | 1800               |

# M-Hilar






Placement for Hilar Cholangiocarcinoma  
Stent-In-Stent Technique  
Uniquely Designed Center Radiopaque Marks  
Second Intervention Possible (PTBD, PDT, etc.)  
Conformability  
Optimal Radial & Returning Force  
Non-Irritating Ends  
Repositionable



## M-Hilar

Coding : B-BONASTENT / B-Biliary / M-M-Hilar  
1-Spiral Type Radiopaque Marks

### 1<sup>st</sup> BONASTENT® M-Hilar




Delivery Device 1800mm Length for T.T.S.

| Order Number | Stent         |                           |                        |                         |                   | Delivery Device |                   |
|--------------|---------------|---------------------------|------------------------|-------------------------|-------------------|-----------------|-------------------|
|              | Diameter (mm) | Proximal Part Length (mm) | Cross Mesh Length (mm) | Distal Part Length (mm) | Total Length (mm) | Diameter (F)    | Usable Length(mm) |
| BBM-100406-1 | 10            | 10                        | 20                     | 10                      | 40                | 7               | 600               |
| BBM-100606-1 | 10            | 15                        | 20                     | 25                      | 60                | 7               | 600               |
| BBM-100806-1 | 10            | 20                        | 20                     | 40                      | 80                | 7               | 600               |
| BBM-101006-1 | 10            | 20                        | 20                     | 60                      | 100               | 7               | 600               |
| BBM-080618-1 | 8             | 20                        | 25                     | 15                      | 60                | 7               | 1800              |
| BBM-080818-1 | 8             | 30                        | 25                     | 25                      | 80                | 7               | 1800              |
| BBM-081018-1 | 8             | 50                        | 25                     | 25                      | 100               | 7               | 1800              |
| BBM-100618-1 | 10            | 20                        | 25                     | 15                      | 60                | 7               | 1800              |
| BBM-100818-1 | 10            | 30                        | 25                     | 25                      | 80                | 7               | 1800              |
| BBM-101018-1 | 10            | 50                        | 25                     | 25                      | 100               | 7               | 1800              |

## M-Hilar


Coding : B-BONASTENT / B-Biliary / M-M-Hilar  
2-Point Type Radiopaque Marks

### 2<sup>nd</sup> BONASTENT® M-Hilar




Delivery Device 1800mm Length for T.T.S.

| Order Number | Stent         |                           |                        |                         |                  | Delivery Device |                   |
|--------------|---------------|---------------------------|------------------------|-------------------------|------------------|-----------------|-------------------|
|              | Diameter (mm) | Proximal Part Length (mm) | Cross Mesh Length (mm) | Distal Part Length (mm) | Total Length(mm) | Diameter (F)    | Usable Length(mm) |
| BBM-100406-2 | 10            | 10                        | 20                     | 10                      | 40               | 7               | 600               |
| BBM-100606-2 | 10            | 15                        | 20                     | 25                      | 60               | 7               | 600               |
| BBM-100806-2 | 10            | 20                        | 20                     | 40                      | 80               | 7               | 600               |
| BBM-101006-2 | 10            | 20                        | 20                     | 60                      | 100              | 7               | 600               |
| BBM-080618-2 | 8             | 20                        | 25                     | 15                      | 60               | 7               | 1800              |
| BBM-080818-2 | 8             | 30                        | 25                     | 25                      | 80               | 7               | 1800              |
| BBM-081018-2 | 8             | 50                        | 25                     | 25                      | 100              | 7               | 1800              |
| BBM-100618-2 | 10            | 20                        | 25                     | 15                      | 60               | 7               | 1800              |
| BBM-100818-2 | 10            | 30                        | 25                     | 25                      | 80               | 7               | 1800              |
| BBM-101018-2 | 10            | 50                        | 25                     | 25                      | 100              | 7               | 1800              |



# Pyloric / Duodenal



- Conformability
- Optimal Radial & Returning Force
- Non-Irritating Ends
- Segmental Compression
- Less Shortening
- Repositionable
- Flare Ends & Partial Covering for Migration Reduction
- Silicone Covering

## Pyloric/Duodenal



Coding : B-BONASTENT / P-Pyloric / Duodenal

Delivery Device 1800 & 2300mm Length for T.T.S.

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BP-200612    | 20            | 60          | 10              | 1200               |
| BP-200812    | 20            | 80          | 10              | 1200               |
| BP-201012    | 20            | 100         | 10              | 1200               |
| BP-200618    | 20            | 60          | 10              | 1800               |
| BP-200818    | 20            | 80          | 10              | 1800               |
| BP-201018    | 20            | 100         | 10              | 1800               |
| BP-200623    | 20            | 60          | 10              | 2300               |
| BP-200823    | 20            | 80          | 10              | 2300               |
| BP-201023    | 20            | 100         | 10              | 2300               |

## Pyloric/ Duodenal Covered





Coding : B-BONASTENT / P-Pyloric / Duodenal / C-Covered

Delivery Device 1800 & 2300mm Length for T.T.S.

| Order Number | Stent         |                     |                   | Delivery Device |                    |
|--------------|---------------|---------------------|-------------------|-----------------|--------------------|
|              | Diameter (mm) | Covered Length (mm) | Total Length (mm) | Diameter (F)    | Usable Length (mm) |
| BPC-200612   | 20            | 30                  | 60                | 10              | 1200               |
| BPC-200812   | 20            | 50                  | 80                | 10              | 1200               |
| BPC-201012   | 20            | 70                  | 100               | 10              | 1200               |
| BPC-200618   | 20            | 30                  | 60                | 10              | 1800               |
| BPC-200818   | 20            | 50                  | 80                | 10              | 1800               |
| BPC-201018   | 20            | 70                  | 100               | 10              | 1800               |
| BPC-200623   | 20            | 30                  | 60                | 10              | 2300               |
| BPC-200823   | 20            | 50                  | 80                | 10              | 2300               |
| BPC-201023   | 20            | 70                  | 100               | 10              | 2300               |

# Gojo-Rectal





Conformability / Optimal Radial & Returning Force  
 Non-Irritating Ends / Segmental Compression  
 Less Shortening / Repositionable / Silicone Covering  
 Flare Ends & Partial Covering for Migration Reduction

## Colo-Rectal



Coding : B-BONASTENT / C-Colorectal

Delivery Device 10F Ø for T.T.S.

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BC-220614    | 22            | 60          | 10              | 1400               |
| BC-220814    | 22            | 80          | 10              | 1400               |
| BC-221014    | 22            | 100         | 10              | 1400               |
| BC-220623    | 22            | 60          | 10              | 2300               |
| BC-220823    | 22            | 80          | 10              | 2300               |
| BC-221023    | 22            | 100         | 10              | 2300               |
| BC-240614    | 24            | 60          | 10              | 1400               |
| BC-240814    | 24            | 80          | 10              | 1400               |
| BC-241014    | 24            | 100         | 10              | 1400               |
| BC-240623    | 24            | 60          | 10              | 2300               |
| BC-240823    | 24            | 80          | 10              | 2300               |
| BC-241023    | 24            | 100         | 10              | 2300               |
| BC-260614    | 26            | 60          | 12              | 1400               |
| BC-260814    | 26            | 80          | 12              | 1400               |
| BC-261014    | 26            | 100         | 12              | 1400               |

## Colo-Rectal Covered





Coding : B-BONASTENT / C-Colorectal / C-Covered

Delivery Device 10F Ø for T.T.S.

| Order Number | Stent         |                     |                   | Delivery Device |                    |
|--------------|---------------|---------------------|-------------------|-----------------|--------------------|
|              | Diameter (mm) | Covered Length (mm) | Total Length (mm) | Diameter (F)    | Usable Length (mm) |
| BCC-220614   | 22            | 30                  | 60                | 10              | 1400               |
| BCC-220814   | 22            | 50                  | 80                | 10              | 1400               |
| BCC-221014   | 22            | 70                  | 100               | 10              | 1400               |
| BCC-220623   | 22            | 30                  | 60                | 10              | 2300               |
| BCC-220823   | 22            | 50                  | 80                | 10              | 2300               |
| BCC-221023   | 22            | 70                  | 100               | 10              | 2300               |
| BCC-240614   | 24            | 30                  | 60                | 12              | 1400               |
| BCC-240814   | 24            | 50                  | 80                | 12              | 1400               |
| BCC-241014   | 24            | 70                  | 100               | 12              | 1400               |
| BCC-260614   | 26            | 30                  | 60                | 12              | 1400               |
| BCC-260814   | 26            | 50                  | 80                | 12              | 1400               |
| BCC-261014   | 26            | 70                  | 100               | 12              | 1400               |

# T raacheal/Bronchial





Conformability  
Performance Formulated Silicone Covering  
Flare Ends for Migration Reduction  
Non-Irritating Ends

## Tracheal/Bronchial



Coding : B-BONASTENT / T-Tracheal / B-Bronchial

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BTB-100309   | 10            | 30          | 8               | 900                |
| BTB-100409   | 10            | 40          | 8               | 900                |
| BTB-100509   | 10            | 50          | 8               | 900                |
| BTB-120309   | 12            | 30          | 9               | 900                |
| BTB-120409   | 12            | 40          | 9               | 900                |
| BTB-120509   | 12            | 50          | 9               | 900                |
| BTB-140309   | 14            | 30          | 10              | 900                |
| BTB-140409   | 14            | 40          | 10              | 900                |
| BTB-140609   | 14            | 60          | 10              | 900                |
| BTB-140809   | 14            | 80          | 10              | 900                |
| BTB-160409   | 16            | 40          | 10              | 900                |
| BTB-160609   | 16            | 60          | 10              | 900                |
| BTB-160809   | 16            | 80          | 10              | 900                |
| BTB-180409   | 18            | 40          | 10              | 900                |
| BTB-180609   | 18            | 60          | 10              | 900                |
| BTB-180809   | 18            | 80          | 10              | 900                |
| BTB-200409   | 20            | 40          | 12              | 900                |
| BTB-200609   | 20            | 60          | 12              | 900                |
| BTB-200809   | 20            | 80          | 12              | 900                |

# Low Stock Items

Inquire for Stock Availability

## Esophageal



Coding : B-BONASTENT / E-Esophageal / R-Reduced Shortening

| Order Number | Stent         |                    |                   | Delivery Device |                    |
|--------------|---------------|--------------------|-------------------|-----------------|--------------------|
|              | Diameter (mm) | Usable Length (mm) | Total Length (mm) | Diameter (mm)   | Usable Length (mm) |
| BE-1806      | 18            | 60                 | 60                | 5               | 700                |
| BE-1816      | 18            | 160                | 160               | 5               | 700                |
| BER-1806     | 18            | 60                 | 60                | 6               | 700                |
| BER-1816     | 18            | 160                | 160               | 6               | 700                |

## Esophageal Shim Anti-Reflux Valve



Coding : B-BONASTENT / E-Esophageal / V-Anti-Reflux Valve

| Order Number | Stent         |                    |                   | Delivery Device |                    |
|--------------|---------------|--------------------|-------------------|-----------------|--------------------|
|              | Diameter (mm) | Usable Length (mm) | Total Length (mm) | Diameter (mm)   | Usable Length (mm) |
| BEV-2009     | 20            | 50                 | 90                | 7               | 700                |
| BEV-2011     | 20            | 70                 | 110               | 7               | 700                |
| BEV-2013     | 20            | 90                 | 130               | 7               | 700                |
| BEV-2015     | 20            | 110                | 150               | 7               | 700                |

## Biliary



Coding : B-BONASTENT / B-Biliary

Delivery Device 1800mm Length for T.T.S.

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BB-080506    | 8             | 50          | 7               | 600                |
| BB-080606    | 8             | 60          | 7               | 600                |
| BB-080706    | 8             | 70          | 7               | 600                |
| BB-080806    | 8             | 80          | 7               | 600                |
| BB-080906    | 8             | 90          | 7               | 600                |
| BB-081006    | 8             | 100         | 7               | 600                |
| BB-080518    | 8             | 50          | 7               | 1800               |
| BB-080618    | 8             | 60          | 7               | 1800               |
| BB-080718    | 8             | 70          | 7               | 1800               |
| BB-080818    | 8             | 80          | 7               | 1800               |
| BB-080918    | 8             | 90          | 7               | 1800               |
| BB-081018    | 8             | 100         | 7               | 1800               |
| BB-100506    | 10            | 50          | 7               | 600                |
| BB-101006    | 10            | 100         | 7               | 600                |
| BB-100518    | 10            | 50          | 7               | 1800               |
| BB-101018    | 10            | 100         | 7               | 1800               |



## Biliary



Coding : B-BONASTENT / B-Biliary / C-Covered

Delivery Device 1800mm Length for T.T.S.

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BBC-100506   | 10            | 50          | 8               | 600                |
| BBC-101006   | 10            | 100         | 8               | 600                |
| BBC-100518   | 10            | 50          | 8               | 1800               |
| BBC-101018   | 10            | 100         | 8               | 1800               |

## Colo-Rectal



Coding : B-BONASTENT / C-Colorectal

Delivery Device 10F Ø for T.T.S.

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BC-221214    | 22            | 120         | 10              | 1400               |
| BC-221414    | 22            | 140         | 10              | 1400               |
| BC-221223    | 22            | 120         | 10              | 2300               |
| BC-221423    | 22            | 140         | 10              | 2300               |
| BC-241214    | 24            | 120         | 10              | 1400               |
| BC-241414    | 24            | 140         | 10              | 1400               |
| BC-241223    | 24            | 120         | 10              | 2300               |
| BC-241423    | 24            | 140         | 10              | 2300               |
| BC-261214    | 26            | 120         | 12              | 1400               |
| BC-261414    | 26            | 140         | 12              | 1400               |

## Colo-Rectal Covered



Coding : B-BONASTENT / C-Colorectal / C-Covered

Delivery Device 10F Ø for T.T.S.

| Order Number | Stent         |                     |                   | Delivery Device |                    |
|--------------|---------------|---------------------|-------------------|-----------------|--------------------|
|              | Diameter (mm) | Covered Length (mm) | Total Length (mm) | Diameter (F)    | Usable Length (mm) |
| BCC-221214   | 22            | 90                  | 120               | 10              | 1400               |
| BCC-221414   | 22            | 110                 | 140               | 10              | 1400               |
| BCC-221223   | 22            | 90                  | 120               | 10              | 2300               |
| BCC-221423   | 22            | 110                 | 140               | 10              | 2300               |
| BCC-241214   | 24            | 90                  | 120               | 12              | 1400               |
| BCC-241414   | 24            | 110                 | 140               | 12              | 1400               |
| BCC-261214   | 26            | 90                  | 120               | 12              | 1400               |
| BCC-261414   | 26            | 110                 | 140               | 12              | 1400               |

# Low Stock Items

Inquire for Stock Availability

## Pyloric/Duodenal



Coding : B-BONASTENT / P-Pyloric / Duodenal

Delivery Device 1800 & 2300mm Length for T.T.S.

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BP-201212    | 20            | 120         | 10              | 1200               |
| BP-201412    | 20            | 140         | 10              | 1200               |
| BP-201218    | 20            | 120         | 10              | 1800               |
| BP-201418    | 20            | 140         | 10              | 1800               |
| BP-201223    | 20            | 120         | 10              | 2300               |
| BP-201423    | 20            | 140         | 10              | 2300               |

## Pyloric/Duodenal Covered



Coding : B-BONASTENT / P-Pyloric / Duodenal / C-Covered

Delivery Device 1800 & 2300mm Length for T.T.S.

| Order Number | Stent         |                     |                   | Delivery Device |                    |
|--------------|---------------|---------------------|-------------------|-----------------|--------------------|
|              | Diameter (mm) | Covered Length (mm) | Total Length (mm) | Diameter (F)    | Usable Length (mm) |
| BPC-201212   | 20            | 90                  | 120               | 10              | 1200               |
| BPC-201412   | 20            | 110                 | 140               | 10              | 1200               |
| BPC-201218   | 20            | 90                  | 120               | 10              | 1800               |
| BPC-201418   | 20            | 110                 | 140               | 10              | 1800               |
| BPC-201223   | 20            | 90                  | 120               | 10              | 2300               |
| BPC-201423   | 20            | 110                 | 140               | 10              | 2300               |

## Tracheal/Bronchial



Coding : B-BONASTENT / T-Tracheal / B-Bronchial

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BTB-100609   | 10            | 60          | 8               | 900                |
| BTB-100709   | 10            | 70          | 8               | 900                |
| BTB-100809   | 10            | 80          | 8               | 900                |
| BTB-120609   | 12            | 60          | 9               | 900                |
| BTB-120709   | 12            | 70          | 9               | 900                |
| BTB-120809   | 12            | 80          | 9               | 900                |
| BTB-140509   | 14            | 50          | 10              | 900                |
| BTB-140709   | 14            | 70          | 10              | 900                |
| BTB-160309   | 16            | 30          | 10              | 900                |
| BTB-160509   | 16            | 50          | 10              | 900                |
| BTB-160709   | 16            | 70          | 10              | 900                |
| BTB-180309   | 18            | 30          | 10              | 900                |
| BTB-180509   | 18            | 50          | 10              | 900                |
| BTB-180709   | 18            | 70          | 10              | 900                |
| BTB-200309   | 20            | 30          | 12              | 900                |
| BTB-200509   | 20            | 50          | 12              | 900                |
| BTB-200709   | 20            | 70          | 12              | 900                |

# Pre-Confirmed Custom Order Items

The items shown below require a Pre-Confirmed Custom Order Basis arrangement. Please, complete the order form to place an order, and inquire with your local representative prior to the issuance of purchase order.

## Esophageal



Coding : B-BONASTENT / E-Esophageal / FG -Flange Type

For the Usage of Upper Esophageal Sphincter (U.E.S.) - **Flange Type**

| Order Number | Stent         |                    |                   | Delivery Device |                    |
|--------------|---------------|--------------------|-------------------|-----------------|--------------------|
|              | Diameter (mm) | Usable Length (mm) | Total Length (mm) | Diameter (mm)   | Usable Length (mm) |
| BEFG-1805    | 18            | 20                 | 50                | 6               | 700                |
| BEFG-1807    | 18            | 40                 | 70                | 6               | 700                |
| BEFG-1809    | 18            | 60                 | 90                | 6               | 700                |
| BEFG-1811    | 18            | 80                 | 110               | 6               | 700                |
| BEFG-1813    | 18            | 100                | 130               | 6               | 700                |
| BEFG-1815    | 18            | 120                | 150               | 6               | 700                |
| BEFG-2005    | 20            | 20                 | 50                | 6               | 700                |
| BEFG-2007    | 20            | 40                 | 70                | 6               | 700                |
| BEFG-2009    | 20            | 60                 | 90                | 6               | 700                |
| BEFG-2011    | 20            | 80                 | 110               | 6               | 700                |
| BEFG-2013    | 20            | 100                | 130               | 6               | 700                |
| BEFG-2015    | 20            | 120                | 150               | 6               | 700                |

## Esophageal



Coding : B-BONASTENT / E-Esophageal / FR -Flare Type

For the Usage of Upper Esophageal Sphincter (U.E.S.) - **Flare Type**

| Order Number | Stent         |                    |                   | Delivery Device |                    |
|--------------|---------------|--------------------|-------------------|-----------------|--------------------|
|              | Diameter (mm) | Usable Length (mm) | Total Length (mm) | Diameter (mm)   | Usable Length (mm) |
| BEFR-1805    | 18            | 20                 | 50                | 6               | 700                |
| BEFR-1807    | 18            | 40                 | 70                | 6               | 700                |
| BEFR-1809    | 18            | 60                 | 90                | 6               | 700                |
| BEFR-1811    | 18            | 80                 | 110               | 6               | 700                |
| BEFR-1813    | 18            | 100                | 130               | 6               | 700                |
| BEFR-1815    | 18            | 120                | 150               | 6               | 700                |
| BEFR-2005    | 20            | 20                 | 50                | 6               | 700                |
| BEFR-2007    | 20            | 40                 | 70                | 6               | 700                |
| BEFR-2009    | 20            | 60                 | 90                | 6               | 700                |
| BEFR-2011    | 20            | 80                 | 110               | 6               | 700                |
| BEFR-2013    | 20            | 100                | 130               | 6               | 700                |
| BEFR-2015    | 20            | 120                | 150               | 6               | 700                |

# Pre-Confirmed Custom Order Items

## Biliary with Retrieval Lasso




Coding : B-BONASTENT / B-Biliary / L-Retrieval Lasso

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BBL-080506   | 8             | 50          | 8               | 600                |
| BBL-080606   | 8             | 60          | 8               | 600                |
| BBL-080706   | 8             | 70          | 8               | 600                |
| BBL-080806   | 8             | 80          | 8               | 600                |
| BBL-080906   | 8             | 90          | 8               | 600                |
| BBL-081006   | 8             | 100         | 8               | 600                |
| BBL-080518   | 8             | 50          | 8               | 1800               |
| BBL-080618   | 8             | 60          | 8               | 1800               |
| BBL-080718   | 8             | 70          | 8               | 1800               |
| BBL-080818   | 8             | 80          | 8               | 1800               |
| BBL-080918   | 8             | 90          | 8               | 1800               |
| BBL-081018   | 8             | 100         | 8               | 1800               |
| BBL-100506   | 10            | 50          | 8               | 600                |
| BBL-100606   | 10            | 60          | 8               | 600                |
| BBL-100706   | 10            | 70          | 8               | 600                |
| BBL-100806   | 10            | 80          | 8               | 600                |
| BBL-100906   | 10            | 90          | 8               | 600                |
| BBL-101006   | 10            | 100         | 8               | 600                |
| BBL-100518   | 10            | 50          | 8               | 1800               |
| BBL-100618   | 10            | 60          | 8               | 1800               |
| BBL-100718   | 10            | 70          | 8               | 1800               |
| BBL-100818   | 10            | 80          | 8               | 1800               |
| BBL-100918   | 10            | 90          | 8               | 1800               |
| BBL-101018   | 10            | 100         | 8               | 1800               |

The items shown below require a Pre-Confirmed Custom Order Basis arrangement. Please, complete the order form to place an order, and inquire with your local representative prior to the issuance of purchase order.

## Biliary Non-Covered in Funnel Shape



Coding : B-BONASTENT / B-Biliary / F-Funnel Shape

| Order Number | Stent         |             | Delivery Device |                    |
|--------------|---------------|-------------|-----------------|--------------------|
|              | Diameter (mm) | Length (mm) | Diameter (F)    | Usable Length (mm) |
| BBF-080506   | 8             | 50          | 7               | 600                |
| BBF-080606   | 8             | 60          | 7               | 600                |
| BBF-080706   | 8             | 70          | 7               | 600                |
| BBF-080806   | 8             | 80          | 7               | 600                |
| BBF-080906   | 8             | 90          | 7               | 600                |
| BBF-081006   | 8             | 100         | 7               | 600                |
| BBF-080518   | 8             | 50          | 7               | 1800               |
| BBF-080618   | 8             | 60          | 7               | 1800               |
| BBF-080718   | 8             | 70          | 7               | 1800               |
| BBF-080818   | 8             | 80          | 7               | 1800               |
| BBF-080918   | 8             | 90          | 7               | 1800               |
| BBF-081018   | 8             | 100         | 7               | 1800               |
| BBF-100506   | 10            | 50          | 7               | 600                |
| BBF-100606   | 10            | 60          | 7               | 600                |
| BBF-100706   | 10            | 70          | 7               | 600                |
| BBF-100806   | 10            | 80          | 7               | 600                |
| BBF-100906   | 10            | 90          | 7               | 600                |
| BBF-101006   | 10            | 100         | 7               | 600                |
| BBF-100518   | 10            | 50          | 7               | 1800               |
| BBF-100618   | 10            | 60          | 7               | 1800               |
| BBF-100718   | 10            | 70          | 7               | 1800               |
| BBF-100818   | 10            | 80          | 7               | 1800               |
| BBF-100918   | 10            | 90          | 7               | 1800               |
| BBF-101018   | 10            | 100         | 7               | 1800               |

# Delivery Devices

## Y – Shape Handle with Braided Tubing

### 10F or Less Delivery Device for

- Biliary Non-Covered & Covered
- Pyloric / Duodenal Non-Covered & Covered
- Colorectal Non-Covered & Covered
- Tracheal / Bronchial

Radiopaque Band



Radiopaque Marking for Distal End of Delivery Device

Yellow Marking Band



Endoscopic View Marking for Proximal End of Captured Stent



## I – Shape Handle with Teflon Tubing

### Larger than 10F Delivery Device for

- Esophageal & Esophageal with Shim Anti-Reflux Valve
- Colorectal Non-Covered & Covered
- Tracheal / Bronchial

Centimeter Marking



Positioning Measurement for Esophageal & A.R.V.





[www.EndoChoice.com](http://www.EndoChoice.com)

 **ENDOCHOICE**  
Tel. +1(770)682 8700  
E-mail. [international@EndoChoice.com](mailto:international@EndoChoice.com)