






PRODUCT NAME	GASTROFLEX <sup>™</sup> UHD	COLOFLEX <sup>™</sup> UHD	CHOLANGIOFLEX <sup>™</sup>	AQ-FLEX <sup>™</sup> 19*	ALVEOFLEX <sup>™</sup>
DISTAL TIP VISUAL					
COMPATIBLE OPERATING CHANNEL	≥ 2.8 mm	≥ 2.8 mm	≥ 1.2 mm	≥ 0.91 mm	≥ 1.9 mm
LENGTH	3 m	4 m	4 m	4m	3 m
MAXIMUM NUMBER OF USES PER PROBE	20	20	10	10	20
USUAL CORRESPONDING PROCEDURE	Eso-Gastro-Duodenoscopy (EGD)	Colonoscopy (compatible with EGD)	Endoscopic Retrograde Cholangio-Pancreatography (ERCP)	Endoscopic UltraSound Fine Needle Aspiration (EUS-FNA)	Bronchoscopy
MAIN PATHOLOGY OF INTEREST	Barrett's Esophagus	- Colonic polyps - Inflammatory Bowel Disease	Indeterminate pancreatico-biliary strictures	Pancreatic cysts	Peripheral Pulmonary Nodules
<b>TECHNICAL SPECIFICATIONS (WITH LSU 488nm)</b>					
FIELD OF VIEW	∅ 240 microns	∅ 240 microns	∅ 325 microns	∅ 325 microns	∅ 600 microns
CONFOCAL DEPTH	55 to 65 microns	55 to 65 microns	40 to 70 microns	40 to 70 microns	0 to 50 microns
RESOLUTION	1 micron	1 micron	3.5 microns	3.5 microns	3.5 microns

\* The AQ-Flex<sup>™</sup> 19 Confocal Miniprobe<sup>™</sup> is CE marked but has not been cleared by the FDA. Product availability cannot be guaranteed in all countries. For further information please contact your local sales representative. Specifications are subject to change without prior notice and without any obligation on the part of the manufacturer.



**Mauna Kea Technologies**  
9, rue d'Enghien  
75010 Paris, France  
Tel +33 1 48 24 03 45  
Fax +33 1 48 24 12 18

[www.maunakeatech.com](http://www.maunakeatech.com)

**Mauna Kea Technologies, Inc.**  
3240 South White Road, Suite 112  
San Jose, CA 95148, USA  
Tel +1 215 478 6288  
Fax +1 215 279 8463