

HALO Cart

Exclusively for HALO Products

Storage space designed to accommodate the HALO Energy Generators, HALO Catheters and accessories.



Compact Design
For working in limited space

Convenient & Efficient
Working Station for the HALO System during the ablation procedure

Custom Design
For HALO Products

HALO Cart

Specifications

A custom procedural cart specifically designed for the storage and easy transport of the HALO System.



Left side view



Back side view



Right side view

Product Information:

HALO Cart: CART-100

Product Features:

- Two shelves accommodate up to two HALO Energy Generators
- Customized trays for HALO³⁶⁰⁺ Catheter products (Holds up to 18 boxes)
- Storage slots for HALO⁹⁰ Ablation Catheters, HALO Guidewires, and HALO Caps
- Tray for HALO footswitch & output cable
- Includes a 1 year warranty from date of purchase

Product Specifications:

- Cart shipped preassembled
- Locking single drawer
- Handles on sides for easy maneuverability
- Light, durable aluminum frame
- Locking wheel casters
- Overall dimensions:
Height: 42.5"
Width: 22.5" (34" w/trays)
Length: 23"
- Net Weight: 100 lbs.

To find out more about the HALO products or to place an order, contact BARRX Medical Toll Free 1.888.66-BARRX (662-2779)

NOTE: Although the complication rate associated with radiofrequency ablation (RFA) using the HALO³⁶⁰ and HALO⁹⁰ ablation catheters has been demonstrated to be low, as with most endoscopic procedures, there are risks associated with RFA, including serious complications. For complete information regarding indications for use, warnings, precautions, adverse events and methods of use, please reference the devices' Instructions for Use.



540 Oakmead Parkway, Sunnyvale, CA 94085
Phone 1.888.662.2779 Fax 408.738.1741
www.barrx.com

© 2010 BARRX Medical, Inc.

L-0102-01 Rev. A