

Accurate, easy-to-use pH monitoring!

Elevate your 24-hour ambulatory pH monitoring to a new level of accuracy with the AccuTrac pHTM monitoring system.

Large LCD screen

4 patient event buttons

Comfortable carrying case



400 - 600% more data recorded and automatically analyzed for every 24 hour session!

Easily adapts to antimony, glass and ISFET pH probes for accurate calibrations and recordings for all patient types.

The AccuTrac pH recorder represents the state of the art in simplicity of use and accuracy in data collection

The compact recorder has a large LCD screen with easy-to-read graphics.

Easy-to-find patient event buttons make symptom reporting and notation of other events, such as supine periods and meal times simple and carefree for the patient.

A lightweight yet durable carrying case protects the recorder during ambulatory use yet still provides access and visibility for all patient interactions during the procedure.

Industry standard USB connections make data upload quick and interference free for maximum data accuracy and efficiency.

*** Also includes buffer solutions for calibrations, calibration stand, test tubes and carrying case.*

