

E.G.Scan

Trans-nasal, disposable system for upper GI screening

E.G.Scan is a system that has been designed to visualise upper GI conditions such as inflammation, Barrett's Oesophagus and ulcers, without the need of patient sedation. The doctor can capture images, has controlled angulations of the probe and produce a report on the findings.

Using their MiroCam capsule endoscopy technology, Intromedic have designed the E.G.Scan system to be small, portable and disposable, giving doctors the freedom to use the system outside of endoscopy or theatres.

The system comprises of 3 parts – A touch screen display monitor, controller and the disposable camera probe.

E.G.Scan – A quick, simple and cost effective tool to diagnose any potential upper GI disorders.

Benefits of the E.G.Scan system

- Disposable
- No sedation required
- Trans-nasal
- Immediate results for patients
- Procedure can be done in clinic
- Report/Analysis
- No obstruction risk

