

## DPU-5000 Series Video Processor

An extraordinary new system providing the best in endoscopic image processing.



- Vivid Full Screen Presentation
- Optional Integrated 15" LCD Display
- Optional Integrated Air Insufflation Control
- Dynamic Image Enhancement Features
- "Pure White" LED Illumination Control
- Complete Stroboscopic Support
- Simple, Portable Image Archiving
- Add-On Camera Support
- Physician Preferred Auto Settings
- Optional Foot Pedal Controls

# A Superior Video Endoscopy System...

## DPU-5000 Series Video Processors

### DPU-5000 Specifications

Image Presentation	Full Screen
Illumination System	LED Control
Formats	NTSC/PAL
Outputs	S-Video (2) RGBS (1) Composite (1) SD Card (1)
Integrated LCD Display	15" Diagonal

### Ordering Information

DPU-5000 Video Processor	07-5000
DPU-5050 Video Processor With integrated LCD Display	07-5050
DPU-5000A Video Processor With integrated Air Insufflation	07-5001
DPU-5050A Video Processor With integrated LCD Display and Air Insufflation	07-5051



### IMAGING FEATURES

- Brightness Controls
- Gain
- Image Enhancement
- Image Freeze with Picture-In-Picture
- Image Orientation Controls
- Full Stroboscopic Support

The DPU-5000 Series of Video Processors present vibrant, sharp CCD based imaging with an extraordinary FULL SCREEN picture. Designed to fit any practice, the DPU-5000 Series has options for any need in today's healthcare environment.

Available with four levels for efficient integration;

- Light-weight, compact processor to work with your existing imaging equipment
- A processor with integrated LCD Display, perfect for the office practice.
- A light-weight, compact processor with integrated air insufflation for esophagoscopy
- A deluxe processor with integrated LCD Display and Air Insufflation; a complete, cost-effective system for any healthcare setting.



### Vision-Sciences, Inc.

40 Ramland Rd. South • Orangeburg, NY 10962

Tel: 845-635-0600 • Fax: 845-635-0601 • E-mail: [info@visionsciences.com](mailto:info@visionsciences.com)

[www.visionsciences.com](http://www.visionsciences.com)

© Vision-Sciences, Inc. Vision Sciences, EndoSheath, and Slide-On are registered trademarks of Vision-Sciences, Inc. Covered by one or more of the following U.S. Patents: 4,825,850; 4,852,551; 4,907,395; 4,947,827; 5,025,778; 5,193,525; 5,271,381; 5,329,887; 5,337,734; 5,443,781; 5,447,148; 5,483,951; 5,518,501; 5,520,607; 5,626,553; 5,667,476; 5,685,822; 5,692,729; 5,702,348; 5,827,177; 5,876,329; 6,174,280; 6,190,330; 6,350,231; 6,461,294; 6,530,881; 6,579,582; 6,733,440; 6,740,030; 6,793,661; 7,025,923; 7,056,284; 7,081,097; Re.34,110. Other U.S. and international patents pending. PL09956 Rev A